


# BE A WINNER IN FIGHT AGAINST COVID-19



## Visual Salocor® SARS-CoV-2 Antigen Test Your Winning Partner in Fight against COVID-19

### Dependable accurate results

- Sensitivity: 96.12%
- Specificity: 99.49%
- LOD: 200 TCID<sub>50</sub>/ml

### Unmatched convenience

- Anterior nasal samples
- Individual buffer vials  
Helping multiple operators at same time
- Visual results  
Easy to interpret
- 15 minutes test time  
Timely results



### Ordering information

Cat.No.	Product	Pack	CE/ EUA Status
COV-07-S	Salocor® SARS-CoV-2 Antigen Rapid Test Cassette	25 T	CE